Biover

Information for BIOVER Disposable Microvascular Clamps  
Informations pour les utilisateurs de clamps microvasculaires BIOVER  
Informationen für die Anwendung mikrovasculärer Klemmen BIOVER  
Informazioni per utenti di clamp microvascolari BIOVER  
Información para los usuarios de las pinzas microvasculares BIOVER  
Informatie voor BIOVER gebruikers van microvaatklemmen  
Information til BIOVER brugere af mikrovasculære klemmer  
Informasjon for BIOVER brukere av mikrokarklemmer  
BIOVER Mikrovasküler Kelepçe Kullanıcıları için Bilgiler  
Πληροφορίες για τους χρήστες των Μικροαγγειακών BIOVER Σφιγκτήρων

## SYMBOL DEFINITIONS:


Do not reuse

Sterilized using irradiation

Do not use if the product sterilization barrier or its packaging is compromised.

CAUTION: Federal (USA) Law restricts this device to sale by, or on the order of, a physician.

Use by date

Batch code

### What are BIOVER disposable microvascular clamps?

BIOVER disposable microvascular clamps are instruments which are used for all microsurgical procedures. The application is to occlude vessels during anastomosis, they may be associated with vessel damage and thrombosis.

### What should the user know about BIOVER disposable microvascular clamps?

These clamps are disposable radiopaque microvascular clamps and can only be used **ONE TIME**.

- The BIOVER disposable microvascular clamps provide an easy identification system for arteries and veins in the surgical field. Both types of disposable microvascular clamps have highly reflective surfaces and are providing reliable deflection under MRI and fluoroscopy.
- BIOVER disposable microvascular clamps deliver a defined and reliable closing pressure. The type of vessel (artery or veins) and the vessel diameter in mm indicates to the surgeon which type of disposable microvascular clamp should be used and which closing pressure is needed to be less traumatic to the vessel.

### How to use BIOVER disposable microvascular clamps?

Before using any type of disposable microvascular clamps, identify the type of vessel and vessel size to make sure that the right clamp with the correct clipping power is used.

- Choose between a single and double BIOVER microvascular clamp.
- Use only products which are in the original and undamaged packaging to ensure that the sterility is guaranteed.
- Check the sterility date on the package and make sure that the date is still valid.

The clamps are gamma sterilized and cannot be resterilized. The manufacturer does not accept any liability for clamps which are resterilized by the user.

- Before closing of the surgical site check that all are disposed.
- After any surgical procedure simple count the empty sterile pouches to account for the precise number of BIOVER disposable microvascular clamps used.

Two types of microvascular clamps are available:

- **Clamps for arteries are yellow**

(Single (1) and Double (2))

- |                    |              |              |
|--------------------|--------------|--------------|
| - Less – 1 mm dia. | Single TKS-1 | Double TKS-2 |
| - 1 mm – 2 mm      | Single TKM-1 | Double TKM-2 |
| - 2 mm – 4 mm      | Single TKL-1 | Double TKL-2 |

- **Clamps for veins are green**

(Single (1) and Double (2))

- |                    |               |               |
|--------------------|---------------|---------------|
| - 0.2 to 1 mm      | Single TKFV-1 | Double TKFV-2 |
| - Less – 1 mm dia. | Single TKSV-1 | Double TKSV-2 |
| - 1 mm – 2 mm      | Single TKMV-1 | Double TKMV-2 |
| - 2 mm             | Single TKLV-1 | Double TKLV-2 |

## WARNINGS

- **Never** use **damaged** clamps and/or **damaged** packaging.
- Use the BIOVER disposable microvascular clamps only **ONE TIME, SINGLE USE**.
- **Do not leave** any BIOVER disposable microvascular clamps in the wound after finishing the surgical procedures.
- **Do not use** the BIOVER vascular clamps for the central blood circulation and not for the central nervous system.

Manufactured by: **BIOVER AG**

Müliweg 2

6052 Hergiswil Switzerland

Tel. **+41 41 630 00 22** • Fax +41 41 630 00 53

e-mail: **biover@solnet.ch**

Distributed by: **Synovis Micro Companies Alliance, Inc**

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