



**THE
DECISION
IS IN
YOUR
HANDS.**

**Keep the power of
nerve information
at your fingertips.**

Information that inspires confidence.

Information that provides insight.

Information that guides decisions.

CHECKPOINT
GUARDIAN™

“Checkpoint directly influences my surgical decision making.”

Scott Kozin, MD | Chief of Staff, Shriners Hospital for Children

WHEN EVERY DECISION IS CRITICAL

In surgery, every decision is critical, and every critical decision must be backed by information. Checkpoint Guardian helps you unlock the information you need to assess and protect nerves during surgery—information that helps you make intraoperative decisions for your patients.

PRECISE ENERGY DELIVERY

Enhanced control and visual confirmation of amplitude and pulse duration output

NERVE SAFE


Biphasic waveform for safe and continuous nerve activation without diminished response

CONFIRMATION OF ENERGY DELIVERY

Two distinct LED indicator lights clearly verify status of stimulation delivery

OBJECTIVE ENERGY OUTPUT

Determine stimulation threshold for motor response throughout the procedure



ORDERING INFORMATION

Product name: Checkpoint Guardian

Description: Surgical Nerve Stimulator/Locator

Instrument Item number: 9095

Unit of measure: Box (4 each) | *RX only*

The Checkpoint Stimulator is a single-use, sterile device intended to provide electrical stimulation of exposed motor nerves or muscle tissue to locate and identify nerves and to test nerve and muscle excitability. Do not use this Stimulator when paralyzing anesthetic agents are in effect, as an absent or inconsistent response to stimulation may result in inaccurate assessment of nerve and muscle function. For a complete list of warnings and precautions regarding the use of the Stimulator please see www.checkpointsurgical.com.

TF: 877.478.9106

O: 216.378.9107

F: 216.378.9116

info@checkpointsurgical.com

Checkpoint Surgical, Inc. 6050 Oak Tree Blvd., Suite 360, Independence, Ohio 44131

CS-9095-MKT-000-A | Dr. Kozin is a paid consultant for Checkpoint Surgical, Inc.



Nerve Care. Empowered.